

TABLE 1 Titration of Hemagglutination (HA) Antigen

	Tube No.						
Reagents (ml)	1	2	3	8	9	10	11 ^a
PBS	0.8	0.5	0.5.....	0.5	0.5	0.5	0.5
Antigen	0.2						
Transfer	0.5→	0.5→	0.5—...→	0.5→	0.5→	0.5→	0.5↓ ^c
0.25% RBC	0.5	0.5	0.5	0.5	0.5	0.5	0.5
Ant. dilution	1:5	1:10	1:20	1:640	1:1280	1:2560	
Results ^b	+	+	+	+	-	-	

^a Tube 11, PBS/RBC control.

^b + = HA; - = no HA (sample titer 1:640).

^c Discard 0.5 ml.